BASTAK
instruments

New Generation

A faint, stylized leaf pattern in the background, consisting of several large, rounded leaf shapes arranged in a symmetrical, branching pattern, rendered in a light gray color against a dark gray background.

*"If you can't measure
it you can't manage it"*


We Produce For
Since

BASTAK

For The World e 1999

Our Founder

Zeki
Demirtařođlu

Our founder Zeki Demirtařođlu started Gazi University Faculty of Agriculture Food Science and Technologies Department after graduating from Ankara Mehmet Rũřtũ Uzel Technical Chemistry High School. After graduating in 1993 he took master education at Food Engineering Microbiology Professorship. Demirtařođlu who became food engineer msc in 1997 worked as quality control manager, business manager and factory manager at various food companies since 1994. He found Bastak Instruments in 1999, started producing flour additives and quality control devices and continuing his works since that day.

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We took our place in developing and globalizing world, we continue to grow.

Bastak instruments which places customer satisfaction on the top of its service pyramid was founded by Zeki Demirtaşoğlu in 1999. Our company started with producing quality control equipments and flour additives, export department was founded in 2003 in line with the objectives, thus took its place in the world market. Bastak instruments which is aiming to comply with the innovations of the globalizing world and be the leader in its field, is making investments for new technologies. Our company, which gives great importance to R&D is aware of necessity of producing standard and high quality products to survive in the world market.

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
Customer satisfaction which is our basic principle is getting more important for us day by day.

Nowadays Bastak instruments is serving to its clients in flour additives, quality control devices for flour and feed, chemical, glass and consumable products for laboratories. Our company is attending national and international fairs as participant or visitor in order to follow developments and introduce innovations. Our company whose basic principle is customer satisfaction is making all investments for new technologies. Our basic values are individual and collective respect, trustworthiness, teamwork, internal and external customer orientation, positive thinking and compromise, self-criticism, understanding and participatory and free expression sharing.

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
GRAIN SAMPLING PROBES

- The device is Bastak branded and its model is 10 000.
- The device, by means of powerful vacuum engine that it has, stores grains such as wheat, barley, paddy, rye and oat as well as heavy products such as corn, bean and chickpea pneumatically (by vacuum) in 2-6 seconds into the sample collection container which is at the laboratory.
- Sample amount can be monitored by means of the sample collection container that has a transparent window on it, and if required, sampling amount can be adjusted.
- Oscillation amount of its powerful and heavy frame is reduced much by means of 4 support sheets of 120 cm.
- Sampling boom of the device can extend up to double hydraulically.
- In this way, the device can scan a field of 8 meters in total as being 4 meters on its right side and 4 meters on its left side.
- The sampling probe can move 120 degrees to the right and 120 degrees to the left, totalling 240degrees.
- The probe of sampler is 2 m.
- The device makes right, left, up, down, forward and back movements hydraulically.
- Compared to electrical systems, the device is much more powerful, fast, silent and long-lived.
- Outside dimensions of the device are 75x75x400 cm.

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NIR GRAIN ANALYZER DA 9000

- The instrument is Bastak brand DA-9000 model.
- It can analysis grains without grinding
- It can analysis wheat within 20 seconds.
- It can make direct analysis on powder products such as flour.
- It can make protein, moisture, ash, gluten and sedimentation analysis.
- There is no need to use any chemical during the analysis.
- The NIR DA-9000 measures between the waves lengths of 908-1694 nm.
- And also it is portable and can be used at field or laboratory.
- A heavy duty rechargeable billet-in battery with the instrument.
- There is no need to purchase different devices for flour and wheat because it can make analyses by replacing only sample container for grain products and powder products.
- It makes data exchange with the computer by means of 2.0 USB input.
- The liquid cooling fan system allows the analysis at the constant temperature.
- It has 10.1 inch sized and 1280x800 pixel resolution touch-screen.
- The outside dimensions of the device are 430x420x410 mm and its net weight is 12 kg.
- The instruments is powered by 220 V or it can be charged and used.

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SAMPLE CLEANER 13500

- BASTAK brand, 13500 model sample cleaner device.
- Through the convertible cylinder sieve inside the device ; sample is sifted and classified . Then sample is collected to the balanced drawers.
- Through the LCD touchable screen, process is able to control, sieve turning speed, feeding gate opening, blowing power, process time and language selection.
- Broken particles, foreign matters such as like powder, feather , good cereals and big foreign matters can be separated easily and classified at wheat, barley, maize, rye, oat etc.
- Before buying raw materials, this device informs about the quality of the products by analysing quality control and impurity.

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Dimensiones: 925x1127x540mm, 85±10kg.

ABSOGRAPH 500

- Comply with the world standards
- 0-200 rpm adjustable speed
- User-friendly, made with care that makes it requires little maintenance
- Temperature control of the kneading Container
- Upon to user use settings changeable by user.
- 21st century most cutting-edge technologies

The instrument can be used for the following applications:

- Grains
- Flour
- Gluten
- Pasta Products
- Dough

Results that can be indicated by the instrument:

- Flour Water absorption
- Viscoelastic of the dough
- Dough constancy and stability
- Dough evolution time
- Dough softening
- Baking behavior
- Gluten index

Related fields:

- Flour Milling
- Bakeries
- Feed mills and Feed factories
- Brewing industry
- Cereal farmers
- Grains Dealers
- Starch factories
- Confectionery /snacks factories

Norm:

- The instrument kneads the flour by using two rotating spiral knives, so the knives faces resistance from the dough, depending on the dough properties, the resistance will be measured by a very accurate seismic analyzer tool, and then a unique coded software invert the incoming data into a graphic chart that makes the data easy to understand and analyze.

Main use:

- The Absograph 500 allows users to understand and analyze
- Water absorption of the flour
- Kneading properties of the flour by following the international standards and norms

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RESISTOGRAPH

- Helps to determine the dough rheology properties by a resistance method.
- Ideal Baking results by determination of the rheology of the dough
- Test your flour improver before using. Results are applicable to real baking process.
- User-friendly and dynamic. Comply with the international standard methods.

The instrument can be used for the following applications:

- Grains
- Flour
- Gluten
- Pasta Products
- Dough

Results that can be indicated by the instrument:

- Flour Water absorption
- Viscoelastic of the dough
- Dough constancy and stability
- Dough evolution time
- Dough softening
- Baking behavior
- Gluten index

Related fields:

- Flour Milling
- Bakeries
- Feed mills and Feed factories
- Brewing industry
- Cereal farmers
- Grains Dealers
- Starch factories

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HECTOLITER

7000

- The device is Bastak branded and its model is 7000.
- It is used for measuring the density of the grains.
- The device is used for measuring the hectoliter in wheat, rice, corn, lentil and similar grains.
- The objective of the hectoliter test in wheat is to estimate the flour productivity.
- When the hectoliter of the wheat increases, the flour productivity increases as well.
- The hectoliter value varies according to the form, density, size and uniformity of the grain. It consists of 6 parts, being measuring container, filling container, transfer container, float, blade and digital balance.
- It is made of draught brass according to world standards. All parts can be separated from each other; they are portable and easy to carry.
- The blade of the device is made of stainless steel.
- The outside dimensions of the device are 180x410x90 mm and its net weight is 3,5 kg.

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NILEMALITRE

7500

- The device is Bastak branded and its model is 7500.
- The objective of the Nilemalitre test in wheat is to estimate the flour productivity.
- When the Nilemalitre of the wheat increases, the flour productivity increases as well.
- The Nilemalitre value varies according to the form, density, size and uniformity of the grain. It consists of 4 parts, measuring container, filling hopper, blade and digital balance.
- It is made of hard brass and Copper coated with chrome, stainless steel and aluminum, according to world standards.
- Portable and easy to carry.
- The blade of the device is made of stainless steel.
- The outside dimensions of the device are 166x132x470 mm and its net weight is 1.75 kg.

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MOISTURE METER 16000

- The device is Bastak branded and its model is 16 000.
- It has digital screen. The product can be selected with the rotary switch
- It measures moisture in 2-5 seconds in flour and wheat and in 12 products in addition to these products.
- Powder and grains products can be measured with the same device.
- The products for which moisture analyse is made: Wheat, rye, barley, oat, rape, corn, beans, durum wheat, sunflower, pea, red wheat, fodder, wheat flour, soybean
- It gives net results since it makes measurement by crushing wheat and other products.
- The outside dimensions of the device are 230 x 150 x 120 mm and its net weight is 2,5 kg.

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HALOGEN MOISTURE ANALYZERS

- Used for quick and precise determination of material humidity upon the basis of weight loss during drying of its small sample (thermo-gravimetric method).
- Halogen radiators ensure effective sample drying.
- Range max 120g, sensitivity 0.001
- Its functioning is based on laboratory scale mechanism.
- It has graphical display and text menu.
- It has drying chart displaying and drying profiles.
- USB, RS232C and PS2 interface
- 4 drying profiles: standard, slow, step and fast
- 20 drying settings memory
- The outside dimensions of the device are 185x290x170 mm and its weight is 4 kg

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4500-s

32

- Bastak branded Roller Mill 4500-s model, double passaged air circulated laboratory type roller semolina mills.
- The device is dumped wheat mill but on will can be used on grinding both dumped and undumped wheat.
- The device consists 2 sections as crushing and grinding in crushing section, the wheat is milled by go through mills and coarse semolina and bran in collector container undersifting part can be separated. On sifting part 850 micron sieve is being used.
- In the grinding section, coarse semolina taking from crushing section can be separated as; flour, fine and coarse semolina. In this section sequentially 160 micron and 280 micronsieves are being used.
- In crushing section there are 3 valces and in grinding section there are 3 valces in total 6 valces are being used.
- In front of the device there are sifting sections. In liso part under sifting section there is 1, in grinding section there is 2 in total 3 collector container is being used.
- In sieving section the lenght of sieve is between 340-360 mm .
- In order to activate the liso part there is 1-0-2 switch button and to activate sieving section there is 0-1 switch button being used. On request it is can run in opposite direction thus the sieve is emptied completely and after switching to other direction can do grinding.
- The device has the capacity of grinding at 100gr/minute.
- The flour efficiency can be up to %40-%70 depending on the wheat quality.
- The outside dimensions of the device are 700x700x900 mm and its weight is 110kgs.
- The device operates at 380V, 50Hz

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ROLLER MILL 4500

- This device is Bastak branded and its model 4500.
- It is a laboratory type double passaged mill.
- It is used to determine the quality of the wheat which will be used for flour production.
- You can get very similar flour like a factory flour
- It is used for both dampened and undampened wheat grinding.
- It is the mixture laboratory mill that enables making necessary modifications and amendments by predetermining the values of the flour to be ground in the factory.
- The device consists of two parts as crushing and liso.
- The wheat milled and divided into three parts as flour, bran(sharp) and semolina by passing through three fluted roller at the crushing section.
- Flour and semolina collected into different drawers by two 160 μ and 800 μ sieves and bran(sharp) taken from the front side of the device.
- Company got informed about the capacity of semolina of wheat right after crushing.
- Semolina which passes through between two flat rollers divided as flour and bran by 160 μ sieve.
- Harmed amylum amount can be controlled by adjusting the distance between soft rollers.
- The device has the capacity of grinding at 600 gr/minute, can be operated 65 %-75% flour efficiency, depending upon the wheat quality.
- Physical, chemical and rheological (water retention and energy graphics) analyses can be done on an obtained flour correctly.
- The outside dimensions of the device are 860x740x650 mm and its net weight is 110 kgs.

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
ROLLER MILL 4000


- The device is Bastak branded and its model is 4000.
- It is a laboratory type roller mill.
- It is used for both dampened and undampened wheat grinding.
- After grinding, flour and bran are taken out from two separate containers, so mixing of them is avoided.
- It is the mixture laboratory mill that enables making necessary modifications and amendments by predetermining the values of the flour to be ground in the factory.
- Depending on the wheat quality, it has a flour productivity of 50–75 %. In the device, grinding speed and amount can be adjusted by means of the feeding adjustment.
- The front part of the device is glass. After grinding, this glass part can be easily removed and the grinding area can be cleaned. In this way, the samples are prevented from mixing into each other.
- The sieving process can be monitored visually by means of the glassed window that is on the section where the sieve unit is located.
- By means of this window, the sieve brush is cleaned after every sample and thus the samples are prevented from mixing into each other.
- Since V belt system is used in the device, no blockage occurs during grinding. The device can operate straightly and reversely.
- The outside dimensions of the device are 400x600x570 mm and its net weight is 42 kg.

Note: Special productions resulting in an extreme quality semolina are performed for macaroni (pasta) factories.

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HAMMER MILL 1900

- Bastak Brand ,1900 Smart Model Hammer Mill.
- Instrument is used to prepare the sample for the Gluten Analysis,Falling Number Analysis and NIR analysis.
- For setting of grinding use air flow rate on it.
- The grinded samples collecting its container.Its motor work 16,800 cycles per minute and can use the diffirent size of sieves.
- Wheat,barley,corn and other grains can grind on it.
- Through the digital touch screen, the user can be operate the rotation and on the touch screen can be monitored on working time.
- After the first sample no need for clean for the second samples.With slow feeding it can grind to the %25 moisture of the grains.
- Sensor on the motor cover when the motor cover is open, dont work for safety.The motor rotational speed : 16,800 rpm.
- Dimensions : 650x570x560 mm . Net weight :51 kg . Grinding capacity : 300 gr in 30-50 seconds.
- The grinding capacity changes depending on the moisture degree.

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HAMMER MILL SMART 1900

- Bastak brand, 1900 Model.
- Instrument is used to prepare the sample for the Gluten Analysis, Falling Number Analysis and NIR analysis.
- For setting of grinding use air flow rate on it.
- The grinded samples collecting its container.
- Its motor work 16,800 cycles per minute and can use the different size of sieves.
- Wheat, barley, corn and other grains can grinding on it .
- After the first sample no need for clean for the second sample. With slow feeding it can grind to the 25% moisture of grains.
- With Air flow regulator can set a flow rate of the samples. Sensor on the motor cover, when motor cover is open dont work for safety. Power is 230 V 50 Hz, the motor rotational speed: 16,800 rpm.
- Dimensions: 640x620x450 mm, net weight: 51 kg, grinding capacity: 300 g in 30-50 seconds.
- The grinding capacity change depending on the moisture degree.

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ENZYME METER 5000

22

- The device is Bastak branded and its model is 5000.
- The device automatically measures enzyme activity in flours and wheats.
- FN measuring mode is used for determining natural alpha amylase enzymes.
- FFN measuring mode is used for determining total (micro-biological + natural) alpha amylase enzymes.
- Changing from FN mode to FFN mode or vice versa can be made by single button.
- The device automatically adjusts the boiling temperature according to the elevation.
- The company name, company address, company phone, fax numbers and webcan be saved in the memory of the device.
- In order to save the results, a printer can be connected to the computer if required and the company information can be seen on every printer output together with the results.
- When the water level reduces, it passes to the stand-by mode for the safety of the operator and a warning message appears on the screen. It has blue LCD screen. The device has 20 function buttons.
- On the screen of the device, date, hour, interior temperature of the device, measuring mode and operation status of the device (running, printing, stop... etc.) can be displayed.
- The outside dimensions of the device are 360x560x440 mm and its net weight is 11 kg.

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ENZYMETER

5100

- The device is Bastak branded and its model is 5100. The device automatically measures enzyme activity in flours and wheats.
- FN measuring mode is used for determining natural alpha amylase enzymes. FFN measuring mode is used for determining total (micro-biological + natural) alpha amylase enzymes.
- The device automatically adjusts the boiling temperature according to the altitude. At the stage of installation company information: company name, company address, company phone, fax numbers and web are saved in the device. The names of the samples to be tested can be entered.
- It automatically calculates the mixing rates in flour and wheat samples. When the moisture values of the samples are entered into the device, it gives the sample amounts to be weighed for testing as (g) grams .
- If the sample amounts are not corrected, it can correct the measuring values obtained as a result of the test according to the moisture.
- It can calculate the liquefaction coefficients of the samples. The device shows both the normal measuring values and corrected values according to the altitude on the printer output.
- The device has blue graphical LCD screen.
- The device has 28 function buttons.
- On the screen of the device, date, hour, ambient temperature, sample names and operation status of the device (running, printing, stop... etc.) can be displayed.
- The outside dimensions of the device are 550x470x190 mm and its net weight is 27 kg.

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6200 - 5200

New Generation

More and More
Instruments are Coming

GLUTEN WASHER 6000

- The device is Bastak branded and its model is 6000.
- It is used to determine the amount of wet gluten in wheat and flour samples by extracting wet gluten.
- The amount of wet gluten of 2 samples can be determined concurrently.
- The device is full automatic.
- If required, in order to determine the amount of wet gluten in special samples, the kneading period of the device can be extended from 20 seconds to 63 seconds and the washing period from 5 minutes to 8,5 minutes.
- The device can also test whole meal.
- On the device, there is "mix" light that is on during the kneading periods and "wash" light that is on during washing periods.
- Moreover, there are "start/stop", "pause" and "reset" buttons.
- The outside dimensions of the device are 320x300x320 mm and its net weight is 23kg.

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GLUTEN WASHER 6100

- The device is Bastak branded and its model is 6100.
- It is used to determine the amount of wet gluten in wheat and flour samples by extracting wet gluten.
- The amount of wet gluten of 2 samples can be determined concurrently.
- If required, in order to determine the amount of wet gluten in special samples, the kneading period and washing periods of the device can be adjusted in the sensitivity of seconds by means of touch buttons on the graphic screen.
- The device can also test whole meal.
- On the device, there is "mix" light that is on during the kneading period, "wash" light that is on during washing periods for the flour samples and "wash meal" light that is on during washing period for the whole meal samples.
- When the moisture values of the samples are entered into the device, it gives the sample amounts to be weighed for testing as (g) grams .
- If the sample amounts are not corrected, it can correct the measuring values obtained as a result of the test according to the moisture.
- The device has blue graphic LCD screen. The device has 27 function buttons.
- On the screen of the device, date, hour, ambient temperature, sample names and the operation status of the device (running, printing, stop... etc.) can be displayed.
- The outside dimensions of the device are 380x370x350mm and its net weight is 22 kg.

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
2100

53

- The device is Bastak branded and its model is 2100.
- Gluten Index is used in determining the quality of wet gluten.
- In this way, it is possible to separate into the groups the samples of wheat and flour according to the qualities of wet gluten.
- Furthermore, the device gives opinion about the sunn pest destruction of the samples of wheat and flour.
- The device is controlled by the micro-processor.
- The device has a blue LCD screen. The device has 6 function buttons.
- On the device's screen, date, hour, environment temperature, test duration, device operating revolution and operation status of the device (ready, testing, stop...etc.) can be displayed.
- The device has to reach to 6000 revolution that is the world standard in 8 seconds and has to operate in 6000 revolution during the remaining 52 seconds period.
- It doesn't operate until the cover is closed due to the security system on its cover, and a warning is displayed on the screen.
- Furthermore, the cover is locked automatically after it starts to operate.
- It automatically stops with the braking system after operating period and gives audible warning.
- The outside dimensions of the device are 200x210x300 mm and its net weight is 8 kg.

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DRY GLUTEN

2500

- The device is Bastak branded and its model is 2500.
- The dry gluten apparatus is used in finding the final and net dry gluten values of the samples of flour and wheat.
- Especially, the values of the dry gluten may be different although the wet gluten values of various wheats and flours having different water absorption are same.
- The device is controlled by micro-processor.
- There are 3 warning lights on the control panel, being power, test on and test off.
- There are thermostat warning light on the device.
- The device automatically adjusts to 150 oC that is the operating temperature.
- When the device reaches to 150 oC, the thermostat warning light turns off.
- The device automatically stops at the end of the testing period of 4 minutes and gives out visual and audible warning.
- The heating surfaces are coated with teflon in order to prevent the samples from adhering. There is a cover lock.
- The outside dimensions of the device are 250x90x200 mm and its net weight is 1,7 kg.

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BREAD VOLUME METER 13300

37

- The Bastak branded model 13300 Volume Testing Device is designed for volumetric testing of oblong and round bakery products.
- The operation of this device is based on the measurement of the volume of millet grains, whose volume is equal to that of the tested bakery product.
- The device is used in laboratories in the bakery industry for quality inspection of bakery products as well as for processing control.
- It works mechanically, thus it requires no electrical source.
- Its ergonomic design enables mobility and it is easy to use.

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Zeleny

3100

- The device is Bastak branded and its model is 3100.
- It is used in determining the bread quality and sunn pest destruction of the samples of wheat and flour.
- It constitutes desired wheat mixture sedimentation and delayed sedimaentation values at the factory.
- The device has graphic LCD screen. There are 5-digit (3 digits for minutes, 2 digits for seconds) 6 separate timers in total on the screen.
- Every timer has a led (light) in order to follow the active timers and it gives light during the operation of the timer.
- The device automatically stops at the end of the period of 5 minutes that is the operating period of the device and gives audible warning.
- The remaining testing period can be seen on the screen during the operation.
- The ambient temperature during testing can be monitored from the screen of the device.
- The device automatically adjusts to 40 revolution/minute that is the world standard whatever the mains voltage is and displays the operating revolution on its screen.
- The outside dimensions of the device are 460x220x350 mm and its net weight is 8 kg.

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LABORATORY SIFTER 8000

8000

- The device is Bastak branded and its model is 8000.
- The device is used for analysing the homogeneity and the particle size of the flour, and the rates of particles to each other.
- Furthermore, it is used in adjusting the distances of the rolls at the flour factories and in controlling the sieving system by supervising the daily production.
- It is with microprocessor control.
- There are power, start and emergency stop buttons and warning lights on the control panel.
- After 5 minutes of testing, the device automatically stops and gives out visual and audible warning.
- When the wooden hoops of the device are replaced with aluminium hoops, in addition to the flour industry, the device can also be used in categorizing the spices by sieving them in the spice industry and in categorizing the solid products such as soil, sand and stone by sieving them in mining and construction industry.
- 7 wooden hoops are supplied with the device.
- 1 collection container is supplied together with these hoops.
- The outside dimensions of the device are 300x580x670 mm and its net weight is 50 kg.

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Muffle

12000

- The device is Bastak branded and its model is 12.000.
- It is used for determining ash amount in the samples of flour and wheat at the flour factories.
- The device is used in the tests intended for burning at high temperature in the industries other than flour industry.
- It has a usable interior volume of 3 dm³.
- It covers very less space at the laboratory. It can analyze 15 samples in one go. Its maximum temperature is 1000 oC.
- Its usage temperature is between 400–950 oC.
- It can operate for 100 hours continuously. Its thermostat adjustment sensitivity is 1 centigrade degree.
- The inside dimensions of the device are 140x110x200 mm and outside dimensions are 320x410x350 mm.
- Its net weight is 20 kg.
- It operates with a power of 1500 W.

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FAT EXTRACTOR Soxholet

It can analyse 6 samples concurrently.

- Temperature can be displayed on the digital screen. It can adjust temperature between 30–200 °C.
- The result can be obtained between 30–60 minutes.
- Since the device execute the extraction process in closed system, its safety is at a higher level than the classical experiment mechanisms. It is generally used in food, feed, detergent, plastic, textile and other industries.
- It provides 60–70 % chemical saving during testing. Six test tubes, extraction container, extraction thimbles, extraction thimble holders are supplied with the device.

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STARCH ANALYSERS Polarimeter

- It is used in starch analyses for raw materials, semi-product and product substances.

- Its measuring range is between 0–180 degrees. Graduation is 1 degree.
- The measuring sensitivity is 0,05 degree.
- It can be used with glass tube of 100 or 200 mm.
- Its light source gives out light at 589 nm. It is generally used in food, grain, feed industries and also in other industries.

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PROTEIN ANALYSER Kjeldahl

- It can burn six samples concurrently. If required, burning units with a capacity of burning twelve and twenty samples can be supplied.
- Sixfold tubespor, exhaust manifolds, digital temperature adjuster and sixfold burning unit are included.
- It takes the caustic in the distillation automatically, as being different from manual model.
- Manual:
 - The device is Raypa branded and its model is DNP-1500, and its code number is 16160015.
 - It can burn six samples concurrently.
 - If required, burning units with a capacity of burning twelve and twenty samples can be supplied.
 - Sixfold tube spor, exhaust manifolds, digital temperature adjuster and sixfoldburning unit are included.

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CELLULOSE ANALYSER Fibre Extractor

- It can analyse raw materials, semi-product and product substances.
- It can analyse 6 samples concurrently.
- It can make 36 analyses in a day.
- It has no chemical loss, the results are absolute.
- Since the device execute the extraction process in closed system, its safety is at a higher level than the classical experiment mechanisms.
- It is used mainly in feed industry and also in food, grain and pharmaceutical industry.
- Six test tubes, extraction case, extraction thimbles, extraction thimble holders are supplied with the device.

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ANALYTICAL Scale 0,0001

- We have different models having a sensitivity of 0,0001
- Without using external grams, the device makes its internal calibration automatically in every start-up.
- Some models have computer and printer connections. In this way, it is possible to save weighing values into an archive.
- There are air windows on the device intended for preventing the effect of the wind. In this way, the result is prevented from being effected in sensitive weighing.
- It is used in ash analyse in flour factories. Moreover, it is used in very sensitive general weighing having a sensitivity of 0,0001.
- All weighing accessories of the device are made of stainless material, and its frame is made of plastic material in order to prevent it from being effected by air conditions and chemical substances.
- The outside dimensions of the device and its net weight change depending on the models.

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SENSITIVE Scale 0,01

- We have different models having a sensitivity of 0,01 gr, from Densi, Shinko brands.
- Some models can be connected to the computer.
- These devices having different sensitivities and capacities can be used in very different enterprises for very different purposes.
- It is used for gluten, sedimentation, index weighing and general weighing at flour factories. The outside dimensions of the device and its net weight change depending on the models.

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DRY AIR STERILIZER OVEN

- Its operating temperature is between the room temperature and 250 C°.
- Its adjustment sensitivity is 1 oC.
- The device can be used for the purposes of sterilization, drying and heating.
- The device is controlled by micro-processor. It has a very useful digital panel for time and temperature adjustments.
- It makes a homogenous heating process.
- It has an electro-acid covered aluminium cell.
- There is a safety thermostat on the device.
- The inside dimensions of the device are 420x320x360 mm and the outside dimensions are 700x470x540 mm. Its net weight is 35 kg.

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WATER DISTILLER

- It has a distilled water capacity of 4/8/12 L/h.
- It has tank capacity 8/16/24L
- It is designed for producing the distilled water required in the laboratories in the quickest and most economical way.
- It uses cooling water of 40 L/h.
- There is a safety sensor on the surface for safety. If it lacks water, it turns off itself automatically.
- The water level in the boiler can be controlled electronically. It has a manometer displaying the pressure of the inlet water.
- All surfaces contacting with water and steam are made of stainless steel.
- It has a warning led for low water pressure, water cut and empty boiler.
- outside dimensions of the device are 810x395x655 mm and its net weight is 28 kg.
- Moreover, it has NS 108 model having a distilled water capacity of 8 L/h and NS 112 models having a distilled water capacity of 12 L/h.

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Purmix C

6-20 g / 50 Kg Flour

- * Improve medium toughness flour
- * Completely getting perfect result at the inefficient kneading time
- * Increase dough processibility
- * Enable soft shell formation

Purmix L

6-20 g / 50 Kg Flour

- * Improve medium toughness flour
- * Increase fermentation tolerance of dough
- * Increase water lift capacity of dough

Purmix M

5-18 g / 50 Kg Flour

- * Enable whiteness of bread clour
- * Increase elasticity of dough
- * Increase bread productivity (ie. increas dough amount being got from unit flour

Purmix P

4-18 g / 50 Kg Flour

- * Increase elasticity of bread shell
- * Enable whiteness of bread colours
- * Increase tolerance of dough
- * Increase process specialty of dough

Purmix R

3-18 g / 50 Kg Flour

- * Increase size and tissue proggress of bread
- * Improve stoma structure of dough and consist of homojen-ic stoma structure
- * Increase proggress specialty of dough
- * Increase gas keeping capacity of dough

Purmix S

3-16 g / 50 Kg Flour

- * Improve weak flour
- * Increase high-tide amount of dough
- * Destroy defaults when dough-kneading
- * Increase water keeping capacity

Purmix T

1-16 g / 50 Kg Flour

- * Improve weak flour
- * Increase swelling amount of dough
- * Increase bread productivity of dough
- * Increase gas keeping capacity of dough

Sünekat AD

8-20 g / 50 Kg Flour

- * Used for improving weak and medium flour
 - * Prevent eurygastrum harm
 - * Increase resistance of dough
- Sünekat AL8-18 g / 50 Kg Flour
- * Used for improving weak and medium flour

- * Prevent eurygastrum harm
- * Increase gas keeping capacity of dough
- * Enable knife-trace of bread

Sünekat AR

4-18 g / 50 Kg Flour

- * Used for improving weak flour
- * Prevent eurygastrum harm
- * Increase fermentation tolerance of dough
- * Enable savings at dough-sour operation

Sünekat AK

6-18 g / 50 Kg Flour

- * Used for improving weak flour
- * Prevent eurygastrum harm
- * Enable extension of kneading time of dough
- * Enable the increase much more bread size
- * Enable knife-trace of bread

Sünekat AS

3-16 g / 50 Kg Flour

- * Used for improving weak flour
- * Prevent eurygastrum harm
- * Put in order stome structure of bread
- * Increase gas keeping capacity of dough
- * Increase energy of dough

Sünekat AT

2-14 g / 50 Kg Flour

- * Used for improving weak and medium flour
- * Prevent eurygastrum harm
- * More stronger gluten structure of dough
- * More stronger protein links of dough
- * Increase gas keeping capacity of dough

Slash 1000

8-18 g / 50 Kg Flour

- * Increase stability of dough
- * Increase fermentation tolerance of dough
- * Prevent eurygastrum harm

Slash 2000

6-15 g / 50 Kg Flour

- * Increase stability of dough
- * Increase fermentation tolerance of dough
- * Increase tolerance against kneading of dough

Slash 3000

4-14 g / 50 Kg Flour

- * Increase stability of dough
- * Increase resistance of dough
- * Increase gas keeping capacity of dough

Slash 4000

2-14 g / 50 Kg Flour

- * Increase stability of dough

- * Increase gas keeping capacity of dough
- * Momentum ripeness of dough

Slash 5000

1-12 g / 50 Kg Flour

- * Momentum ripeness of dough
- * Increase gas keeping capacity of dough
- * Increase bread size and prolong shelf-life

Slash 6000

1-10 g / 50 Kg Flour

- * Momentum ripeness of dough
- * Being homogen stoma structure of bread
- * Increase gas keeping capacity of dough
- * Prolong shelf-life of bread

Arkat E10

0,5-18 g / 50 Kg Flour

- * Increase kneading tolerance of dough
- * Increase elasticity of dough
- * Increase process of dough

Arkat E25

0,5-15 g / 50 Kg Flour

- * Increase kneading tolerance of dough
- * Increase elasticity of dough
- * Increase size of bread

Arkat E40

0,5-13 g / 50 Kg Flour

- * Used for medium strength flour
- * Increase energy value of dough
- * Decrease shell thickness of bread
- * Enable kneading time shortening

Arkat E50

0,5-10 g / 50 Kg Flour

- * Used for medium strength flour
- * Prevent little eurygastrum harm
- * Enable kneading time shortening
- * Minimize process defaults at the oven

Arkat E60

0,5-8 g / 50 Kg Flour

- * Used for medium strength flour
- * Consist of small homogen gas rooms in bread
- * Enable discharge of bread's inside at the oven step
- * Gain soft structure to bread

Arkat E100

0,5-6 g / 50 Kg Flour

- * Increase size and stoma improvement of bread
- * Improve stoma structure of dough and become homogen stoma structure
- * Increase process specialty of dough
- * Increase gas keeping capacity of dough

Armix 5000

5-20 g / 50 Kg Flour

- * Increase elasticity of dough
- * Increase swelling amount
- * Enable loosness of dough
- * Prevent dough drying

Armix 6000

5-18 g / 50 Kg Flour

- * Increase elasticity of dough
- * Increase swelling amount
- * Prevent torn and broken pieces
- * Prevent cracking and being poured out at the bread-shell

Armix 10000

3-17 g / 50 Kg Flour

- * Increase operation of strong flour
- * Decrease resistance of dough
- * Increase gas amount formation in bread
- * Enable kneading time shortening
- * Enable brawn-shell with bright lively colours

Armix 60000

2-16 g / 50 Kg Flour

- * Increase operation of strong flour
- * Increase gas amount of dough
- * Decrease shell thickness of bread
- * Balance inefficient enzyme activity
- * Improve lengthenning speciality.

Armix 80000

2-14 g / 50 Kg Flour

- * Increase elasticity of dough
- * Increase swelling amount
- * Decrease resistance of dough
- * Prevent the getting stale and braking at bread

Armix 105000

1-10 g / 50 Kg Flour

- * Increase elasticity of dough
- * Increase swelling amount
- * Enable brawn-shell with bright lively colours
- * Prevent the getting stale and braking

Armix 120000

1-18 g / 50 Kg Flour

- * Increase very strong process of flour
- * Enable brawn-shell with bright lively colours
- * Decrease resistance of dough
- * Prevent the getting stale and braking
- * Increase specialty being given easily to dough

BASTAK
instruments

**WE PRODUCE
FOR THE
WORLD**

Office / Ofis

İvedik OSB 1354. Cad. 1387. Sk. No: 19
Yenimahalle/Ankara/TURKEY
Tel: +90 312 395 6787 (pbx) Fax: +90 312 395 6788
E-mail: export@bastak.com.tr

Factory / Fabrika

Dağyaka Mah. 2008. Sok. No: 1
Kahramankazan/Ankara/TURKEY
Tel: +90 312 811 2811 (pbx) Fax: +90 312 811 2812
E-mail: export@bastak.com.tr

www.bastak.com